

# Microbiology Focus

Volume 5.4, 2013

**Fluka**  
Analytical

## *Clostridium perfringens*: A Pathogen and Fecal Indicator



*Clostridium perfringens* is an indicator for fecal contamination and is also a pathogen. Water is one possible source for further contamination as the spores are quite resistant to heating and chlorination. A scanning electron microscope image of vegetative and spore forming cells from *C. perfringens*.

(source *C. perfringens* image: Copyright Dennis Kunkel Microscopy, Inc.)

Fung Double Tube Method .....	2
<i>Clostridium perfringens</i> Chromogenic Media .....	3
Bifidobacteria .....	6

# Fung Double Tube Method

Daniel Y. C. Fung, Food Science, Kansas State University; Roger Fujioka, Professor, University of Hawaii; Josep Yuste and Marta Capellas, Veterinary College, Universitat Autònoma de Barcelona — [dfung@k-state.edu](mailto:dfung@k-state.edu)

## Four-hour rapid evaluation of fecal contamination of recreational water in seaside beaches by the Fung Double Tube Method.

Recreational waters such as those found at seaside beaches, rivers, lakes and ponds, are key areas for recreational use such as swimming, bathing, splashing, water games, etc. However, many of these waters may be contaminated with discharge, sometimes with fecal material from towns, cities, and small and large municipalities. These waters may contain fecal bacteria and other microbes which may pose a public health risk for swimmers, bathers and the general public who use these waters. Thus, it is important to monitor these waters regularly for the presence or absence of fecal contamination to protect the safety of the swimming public.

Conventional detection of fecal contamination in water takes one to two days. By the time the results are available, the quality of the water being tested may have changed substantially, making the data of very limited value.

In the past two years (2011–2013), the Water Research Group, a division of the Department of Public Health in Honolulu, Hawaii, under the direction of Dr. Roger Fujioka, has been testing newer and faster methods to ascertain potential fecal contamination of beach water. Using the unique Fung Double Tube Method, which can generate four-hour data from the time the water is applied to the system to the time results are read, allows a positive or negative result for fecal contamination. Currently, this is the fastest known method in the world for detection of fecal bacteria in water.

The heart of the test is as follows. *Clostridium perfringens*, an anaerobic bacterium which produces spores, is ubiquitous in fecal material (animal and human). Thus, the presence of *C. perfringens* in seawater would indicate a high probability of fecal contamination. The conventional method to detect the presence of *C. perfringens* in water takes 24 to 48 hours, making the results of very little value to warn the public of the potential danger of the water. Daniel Y. C. Fung, Ph.D., Professor of Food Science at Kansas State University, developed a Double Tube Method which can detect and enumerate *C. perfringens* in water in only four hours, making this the fastest test for fecal indicator organisms in the world. The reason for this success is that *C. perfringens* has the remarkably fast generation time of only 7.1 minutes at 41 °C, making it the fastest growing bacterium known to microbiologists at this time.

The group in Honolulu has been testing the Fung Double Tube System for about two years, and has found the method to be clean and very easy-to-use. It can indeed obtain four-hour results so that decisions can be made in four hours to open or close the beach for recreational use.

## Did you know...

**The production of *Clostridium perfringens* toxins starts under certain conditions?**

Only one serotype out of five produces toxins, and the toxin production doesn't start until a concentration of  $10^8$  vegetative germs per gram of food or beverage is reached.



Figure 1: Keeping a pot simmering between 5–65 °C can provide good conditions for the growth of *C. perfringens*.





Figure 2: The Fung Double Tube System (for cultivation of anaerobes aerobically)



**Figure 3:** Each large tube contained 10 mL of the water sample selected for experimentation added into 10 mL of double strength (2X) melted sterile SFP agar. After rotating the larger tube with the sample and melted agar for 10 seconds to mix, a sterile inner tube is carefully inserted into the larger tube, thereby squeezing the liquid sample and SFP agar between the space of the two tubes, forming a thin film of sample and agar between the two tubes. The cap being tightly screwed, as well as a very tight, thin agar, makes the system very anaerobic, ideal for growth of *C. perfringens*. After about four hours of incubation at 41 °C, each single cell is visible as a black colony. *C. perfringens* is the fastest growing bacterium known; the generation time, under ideal conditions, is only 7.1 minutes.

The Fung Double Tube System involves two tubes—the smaller inside tube and a larger outer tube. Into the larger tube, 10 mL of liquefied Shahadi Ferguson Perfringens (SFP) agar (specially formulated to detect *C. perfringens*) was first placed and then sterilized. The medium is kept at 42 °C so that it will be in a liquid form. When a water sample is to be tested, 10 mL of the water (seawater, river water, etc.) will be carefully introduced into the large tube with the SFP agar. Then, the final step is to insert a smaller tube into the larger tube which contains both the SFP agar and the water sample. As the smaller tube is inserted into the large tube, the medium and sample will become sandwiched between the two tubes and form a thin film in the cavity. The material will be pushed up to the tip of the two tubes, then a screw cap will be applied, making this a very anaerobic system. In four hours, if there are *C. perfringens* in the water sample, tiny black colonies will form. The colonies will grow in size as time progresses, but at four hours, a researcher can start to ascertain if the water does or does not have *C. perfringens* and how many per mL of water.

If there are no black colonies found in the water, then it can be determined that the tested water is free of fecal material and safe to use for recreational purposes. If the number of colonies is 1–10, then the contamination level would be small and may be safe for swimming with caution. However, if the number is 10–100 colonies, the water is not safe for swimming, and the beach should be closed.

## Rapid Detection of *Clostridium perfringens* by a New Chromogenic Media

Mohammad Manafi, Institute for Hygiene and Applied Immunology, Medical University of Vienna;  
Jvo Siegrist, Product Manager- Microbiology — [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)

### Selective chromogenic media for isolation and enumeration of *Clostridium perfringens* in water samples using membrane filtration.

CP *ChromoSelect* Agar is a selective chromogenic media for isolation and enumeration of *Clostridium perfringens* in water samples using membrane filtration. This new agar is more reliable and easier to handle than m-CP and TSC agars. The color does not diffuse in the agar and confirmation is not required since the green coloration is specific for *C. perfringens*. *C. perfringens* is an anaerobic, Gram-positive, spore-forming rod-shaped bacterium. It is widespread in the environment and is also found in the digestive systems of humans and domestic and feral animals. Perfringens poisoning, usually the result of ingesting under-cooked food, especially meat, is one of the most commonly reported foodborne illnesses. Early detection of *Clostridium* in food and water is important to control outbreaks. *C. perfringens* produces an extensive range of invasins and exotoxins. The enterotoxins cause the undesirable, mostly meat-associated food poisoning, and wound and surgical infections that lead to gas gangrene. *C. perfringens* plays a subsidiary role in water examination.<sup>2</sup> Clostridia are spore builders and are resistant to heating, chlorination and other stress factors. In contrast to vegetative cells like coliforms

### Did you know...

#### *Clostridium perfringens* is a special anaerobe?

It is a strictly anaerobic bacterium but is able to survive when exposed to oxygen for short periods of time. A complex adaptive response to reactive oxygen species was observed but not completely understood.




**Figure 1:** A photomicrograph of a Gram-stained culture specimen from a patient with gas gangrene, showing numerous *Clostridium perfringens*. (Source: CDC 1979).




(*E. coli*, enterococci), which are less resistant, *C. perfringens* has the advantage of surviving longer.<sup>3</sup> Therefore, while fecal contamination is detected mostly by coliforms as an indicator, which could disappear after a processing step, *C. perfringens* remains present. The organism is not a hazard in water; rather, it is problematic when the water comes in contact with food. In consideration of the aforementioned facts, it is obvious that detection and identification of *C. perfringens* is an important step toward the control and eradication of this potent pathogen. Some characteristic enzymes of *C. perfringens* are: hemolysins ( $\beta$ -hemolysis), lecithinase, extracellular proteases, lipases (phospholipase-C), collagenase, hyaluronidase, saccharolytic, and enzymes to reduce sulphite to sulphide. These enzymes are also used as detection and differentiation targets. It is also notable that *C. perfringens* is a non-motile bacterium, and it is the most important of the sulphite reducing clostridia. *C. perfringens* normally grows at 44 °C, whereas some other clostridia are inhibited at this temperature. This property is used in ISO methods to give the medium more selectivity.<sup>4</sup>

Early detection of *Clostridium* in food is important to control outbreaks. To facilitate detection, we have introduced a new chromogenic media, CP *ChromoSelect* Agar, for enumeration and differentiation of *Clostridium* sp., in particular *Clostridium perfringens*, in aqueous samples. In the present study, three media types (mCP, TSCF and CP *ChromoSelect* Agar) were evaluated for recovery of *C. perfringens* in different surface water samples. Using a membrane filtration technique on 139 water samples, 131 samples (94.2%) were found to be positive for *C. perfringens* in at least one of the culture media. Green colored colonies on CP *ChromoSelect* Agar (Figure 2) were counted as presumptive *C. perfringens* isolates.



**Figure 2:** Drinking water sample cultured on CP *ChromoSelect* Agar. *C. perfringens* appears as distinct green colonies.




**Figure 3:** *C. perfringens* ATCC 10873 on CP *ChromoSelect* Agar (left) and TSC agar (right) (Note the false negatives on the TSC agar).

For detection of *C. perfringens*, m-CP and TSC agar have been recommended.<sup>2,5</sup> However, there are problems associated with each of these media. CP *ChromoSelect* Agar is more reliable and easier to handle than m-CP and TSC agars. The color does not diffuse in the agar and confirmation is not required since the green coloration is specific for *C. perfringens*. In addition, the recovery of *C. perfringens* was rejected by ISO in favor of methods based on TSC agar.<sup>3,6</sup> CP *ChromoSelect* Agar also eliminates the excessive and variable blackening of the peripheral colonies encountered with TSC agar, which makes colony counting at lower dilutions difficult, and leads to false positives. It is also more reliable at high bacteria counts, where the TSC agar can produce false negatives because of interference with the other enzymatic mechanisms from acid production and oxygen contact (Figure 3). TSC detects all sulphite-reducing clostridia, and not only *C. perfringens*.

Out of 483 green colonies on CP *ChromoSelect* Agar, 96.3% (465 strains, indole negative) were identified as *C. perfringens*, 15 strains (3.1%) were indole positive and were identified as *C. sordelli*, *C. bifermentans* or *C. tetani*. Only 3 strains (0.6%) gave false positive results and were identified as *C. fallax*, *C. botulinum*, and *C. tertium* (Table 1). Variance analysis of the obtained data showed no statistically significant differences in the counts obtained between media used in this work (Figure 5). In general, the identification of typical and atypical colonies isolated from all media demonstrated that CP *ChromoSelect* Agar was the most useful medium for *C. perfringens* recovery in water samples.

CP *ChromoSelect* Agar avoids the disadvantages of m-CP agar, such as the presence of ammonia that prevents subculturing the *C. perfringens* colonies, the too-selective nature of m-CP agar and the evanescence of the red color of colonies after the addition of ammonia, which makes further confirmation impossible (Figure 4).




Strains	Indole Reaction	n
<i>C. perfringens</i>	-	465 (96.3%)
<i>C. tertium</i>	-	1 ( 0.2%)
<i>C. botulinum</i>	-	1 ( 0.2%)
<i>C. fallax</i>	-	1 ( 0.2%)
<i>C. bifermentans</i>	+	2 ( 0.4%)
<i>C. sordelli</i>	+	12 ( 2.5%)
<i>C. tetani</i>	+	1 ( 0.2%)

**Table 1:** All green colonies isolated from CP *ChromoSelect* Agar and identified with API system (n=483)


Besides its advantages over m-CP and TSC agars, CP *ChromoSelect* Agar is an ideal growth media. It contains only vegetable peptones and, together with yeast extract, it is an excellent source of nitrogen, carbon, amino acids and vitamin B complex. Sucrose acts as the fermentable carbohydrate, and reducing agents lower the redox potential of the media. Diverse salts provide the required ions for enzymatic reactions. Buffering agents stabilize the pH within the ideal growth range. Inhibitors Dcycloserine and polymyxin B give the medium its selectivity, while further selectivity is achieved by incubation under anaerobic conditions at 44 °C. Various promoters and substrates protect injured cells to improve recovery rate and enhance growth. The chromogenic enzyme substrates in the CP *ChromoSelect* Agar provide the differentiation, for *C. perfringens* in particular (Table 2). A negative indole reaction (Kovac's Reagent) is confirmatory for *C. perfringens*.

Organisms (ATCC)	Growth	Colony Appearance
<i>Clostridium perfringens</i> (13124)	+++	Green
<i>Clostridium bifermentans</i> (638)	+++*	Dark blue with violet halo
<i>Clostridium sporogenes</i> (8534)	-	-
<i>Clostridium sordelli</i> (9714)	++	Dark green with halo (change to red with Kovac's Reagent)
<i>Enterococcus faecalis</i> (29212)	++	Violet
<i>Escherichia coli</i> (25922)	-	-
<i>Pseudomonas aeruginosa</i> (27853)	-	Colorless
<i>Staphylococcus aureus</i> (25923)	-	-
<i>Bacillus subtilis</i> (6051)	-	-
<i>Salmonella typhimurium</i> (DSM 554)	++	Violet

**Table 2:** *Clostridium* sp. cultural characteristics in CP *ChromoSelect* Agar.  
\*Growth at 40 °C, but no growth at 44 °C



**Figure 4:** *C. perfringens* on m-CP Agar



**Figure 5:** Comparison between TSCF agar and the m-CP and CPC (CP *ChromoSelect* Agar) media for enumerating strains of *C. perfringens* in water samples.

- CP *ChromoSelect* Agar was the most useful medium for *C. perfringens* recovery in water samples.
- CP *ChromoSelect* Agar is more reliable and easier to handle than m-CP and TSC agars. The color does not diffuse in the agar and confirmation is not required since the green coloration is specific for *C. perfringens*.

#### References

1. Mammad Manafia, Kerstin Waldherra, Michael Kundib, Evaluation of CP *ChromoSelect* Agar for the enumeration of *Clostridium perfringens* from water, International Journal of Food Microbiology (2013).
2. ISO 6461 Part 2. Water Quality – detection and enumeration of *Clostridium perfringens* (Revision of ISO 6461-2: 1986).
3. National Standard Method, Enumeration of *Clostridium perfringens* by Membrane Filtration, issue no. 3.1, reference no W5i3.1 (2005).
4. Wells, C. L.; Wilkins, T. D., Clostridia: Sporeforming Anaerobic Bacilli, In Baron's Medical Microbiology, 4th ed., S. Baron Editor, Univ. of Texas Medical Branch, 1996.
5. Council Directive 98/83/EC, Relating to the Quality of Water Intended for Human Consumption, Official Journal of the European Communities, p.L330-32-L330/53 (1998).
6. Sartory, D. P.; Field, M.; Curbishley, S.M.; Pritchard, A.M. Evaluation of Two Media for the Membrane Filtration Enumeration of *Clostridium perfringens* from Water. Lett. Appl. Microbiol. 1998, 27, 323–327.

Sigma-Aldrich's range of media and supplements for *Clostridium perfringens* detection

Media	Description	Cat. No.
<b>Nonselective Media</b>		
AEA Sporulation Broth (Base), modified	For early sporulation of <i>C. perfringens</i> from foods.	17170
<b>Nonselective Differential Media</b>		
Gelatin Iron Agar	Used for detecting gelatin liquefaction and hydrogen sulphide production.	G0289
Lactose Gelatin Medium (Base)*	For the detection of lactose and gelatine metabolizing microorganisms ( <i>C. perfringens</i> ).	61348
Motility Nitrate Medium	Selective medium for motile nitrate-utilizing microorganisms ( <i>Ps. aeruginosa</i> , <i>C. perfringens</i> ).	14305
<b>Selective Differential Media</b>		
CP ChromoSelect Agar	A chromogenic agar for enumeration, detection and differentiation of <i>Clostridium</i> sp., in particular <i>C. perfringens</i> , in aqueous samples.	12398
m-CP Agar Base	Recommended by the Directive of the Council of the European Union 98/83/EC for isolation and enumeration of <i>C. perfringens</i> from water samples using membrane filtration technique.	75605
Perfringens Agar Base	Used for the presumptive identification and enumeration of <i>C. perfringens</i> from food.	39727
SPS Agar	Medium for detection, isolation and enumeration of <i>C. perfringens</i> and <i>C. botulinum</i> in food acc. to Angelotti et al. (1962).	85627
SPS Agar, modified	For the selective isolation and enumeration of <i>C. perfringens</i> from foods.	17231
TSC Agar* (Base for TSC Fluorescence Agar, SFP Agar, etc.)	For the isolation and enumeration of vegetative forms as well as spores from <i>C. perfringens</i> in food, clinical specimens and other material.	93745
TSN Agar*	Highly selective medium for the detection and enumeration of <i>C. perfringens</i> in food and other material.	93735

\*not available in U.S.A.

Supplements	Cat. No.
Egg Yolk Emulsion	17148
m-CP Selective Supplement I	51962
m-CP Selective Supplement II	82265
Perfringens S.F.P. Selective Supplement	53436
Perfringens T.S.C. Supplement	P9352
TSC Agar Supplement	80548

***Clostridium perfringens* as Certified Reference Material**

Evaluate the performance of your test with a certified microorganism standard.

Vitroids™ Test Strain	Origin	Strain No.	CFU	Cat. No.
<i>Clostridium perfringens</i>	NCTC	10240	30	RQC02351
			200	RQC02355
			500	RQC20106

## Bifidobacteria in the Dairy Industry

Jvo Siegrist, Product Manager- Microbiology — [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)



### Bifidobacteria and the quality control of the dairy products they produce.

Bifidobacteria like *Bifidobacterium longum*, *infantis* and *brevi* are used for manufacturing dairy products. Bifidobacteria are Gram-positive, non-motile, rod-shaped and often branched obligate anaerobic bacteria. As probiotics, they have a positive affect on the immune system, control the intestinal pH and the gut microflora. Bifidobacteria produce bacteriocins and bacteriocin-like inhibitory compounds which hinder the growth of undesirable bacteria in the intestine. *Bifidobacterium* species were discovered by Tissier in 1889, and typically came to be associated with feces, especially those of breast-fed infants.<sup>2</sup>

*B. longum* is the best characterized species in the genus *Bifidobacterium*. It is able to utilize a broad range of substrates for energy, such as plant polymers, glycoproteins and glycoconjugates, as well as having specialized proteins for the catabolism of oligosaccharides like galactooligosaccharide

(GOS) and fructooligosaccharide (FOS). GOS and FOS are prebiotics and are also often used as dietary supplements. A so-called bifidogenic growth stimulator (BGS) produced by *Propionibacterium freudenreichii* was observed in Swiss cheese, which improved the health of patients with ulcerative colitis. This was likely a prebiotic effect as well.<sup>3</sup>

Bifidobacteria also have a unique hexose metabolism called the bifid shunt. The key enzyme, fructose-6-phosphate phosphoketolase, is not found in any other Gram-positive intestinal bacteria and therefore provides an ideal target for a diagnostic test. The optimum growth temperature is 37–41 °C and they produce both lactic (L-isomer) and acetic acids from lactose. Sodium propionate is used to promote the growth of bifidobacteria. Cysteine helps in creating reduced conditions required for the growth of anaerobic bifidobacteria. *Bifidobacterium* species are known to be acid and bile tolerant. Other known antibiotics used for selective media are mupirocin, kanamycin B, polymyxin B, neomycin, paromomycin, gentamicin and nalidixic acid.




**Figure 1:** Scanning electron microscope image of *Bifidobacterium*  
(Source: Bao Qiuhua Inner Mongolia Agricultural University)

In an adult's intestine, only 3–6% of the fecal flora is composed of bifidobacteria, while in breast-fed infants, bifidobacteria can be up to 90%. With increasing age, the number of bifidobacteria decrease. Human milk contains many substances that stimulate the growth of bifidobacteria *in vitro*. It was observed that babies and adults with lower numbers of bifidobacteria have a higher risk for diarrhea and allergies. For this reason, *Bifidobacterium* are added as a probiotic supplement to infant formulas, drinks, yogurts and many other products.

### Bifidobacteria Selective Media (BSM) with color reaction

Because of the wide use of bifidobacteria and a customer's request, Sigma-Aldrich has developed a Bifidobacteria Selective Media (BSM), available as an agar or a broth, for daily quality control. This medium allows for fast and easy quality control of yogurt made with bifidobacteria and can be used to control the count of bifidus bacteria. *Bifidobacterium* grows very well on this medium, while *Lactobacillus* and *Streptococcus* strains are inhibited. *Bifidobacterium* colonies grow within 24–48 hours (occasionally up to three days because of the highly selective conditions). The *Bifidobacterium* colonies are purple-brown and therefore easy to differentiate from other organisms.



**Figure 2:** Yogurt sample cultured on BSM Agar. Bifidobacteria appear as purple-brown colonies.

In a Swiss governmental evaluation study for the enumeration of bifidobacteria in sour milk products, the traditional method was compared to Wilkins-Chalgren Agar with 100 mg/L mupirocin and BSM Agar. The traditional method gave significantly different results, while Wilkins-Chalgren Agar and BSM Agar showed similar results without any significant differences. The study remarked: "On the BSM Agar, the bifidobacteria forms purple-brown colonies which made the enumeration easy".<sup>1</sup>

Organisms (ATCC)	Growth	Colony Appearance
<i>Bifidobacterium longum</i> (15707)	+++	red-brown (maroon)
<i>Bifidobacterium infantis</i> (15697)	+++	red-brown (maroon)
<i>Streptococcus thermophilus</i> (14486)	–	–
<i>Lactobacillus acidophilus</i> (314)	–	–
<i>Lactobacillus bulgaricus</i> (11842)	–	–

**Table 1:** *Bifidobacterium* sp. cultural characteristics on BSM Agar

### TOS-propionate agar

The TOS-propionate agar medium is recommended by ISO/IDF for enumeration of presumptive bifidobacteria by colony count technique from milk products with mixed flora. The medium is selective and other lactic acid bacteria are inhibited. In addition, the medium contains galactooligosaccharides, which is one of the most excellent bifidobacteria growth promoting substances. Lithium mupirocin is the antimicrobial supplement that inhibits the growth of most lactic acid bacteria commonly used in fermented and non-fermented dairy products.<sup>4,6</sup>

Name	Cat. No.
Bifido Selective Supplement B	90577-5VL
BSM Agar	88517-500G-F
BSM Broth	90273-500G-F
BSM Supplement	83055-5G-F
Lithium mupirocin Supplement	69732-10X5ML 69732-10X25ML
TOS-propionate Agar	43314-500G
Wilkins Chalgren Anaerobic Agar	W1761-500G

**Table 2:** Selective medium and supplements for *Bifidobacterium*

### References

1. Bundesamt für Gesundheit: Schweizerisches Lebensmittelbuch (SLMB): Kapitel 56, «Mikrobiologie» Neuauflage 2000, Stand 2004.
2. Tissier H., Recherches sur La Flore Intestinale. Normale et Pathologique du Nourrisson., Georges Carre et C. Naud, Eds., 3 Paris (1900).
3. Suzuki A., et al., Bifidogenic Growth Stimulator for the Treatment of Active Ulcerative Colitis: A Pilot Study, *Nutrition*, 22(1), 76–81 (2006).
4. ISO/DIS 29981 IDF 220, Milk products- Enumeration of presumptive bifidobacteria- colony count technique at 37 °C, 2010.
5. Rada V. and Koc J.: The Use of Mupirocin for Selective Enumeration of Bifidobacteria in Fermented Milk Products. *Milchwissenschaft*. 55, 65–67 (2000).
6. Thitaram, S.N., Siragusa G.R., Hinton A., *Letters in Applied Microbiology*, Vol 41, 355-360, Bifidobacterium Selective Isolation and Enumeration from Chicken Ceca by an Oligosaccharide- antibiotic Selective Agar Medium (2005).



## Sigma-Aldrich® Worldwide Offices

### Argentina

Free Tel: 1800 888 7446  
Tel: (+54) 11 4556 1472  
Fax: (+54) 11 4552 1698

### Australia

Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Austria

Tel: (+43) 1 605 81 10  
Fax: (+43) 1 605 81 20

### Belgium

Tel: (+32) 3 899 13 01  
Fax: (+32) 3 899 13 11

### Brazil

Free Tel: 0800 701 7425  
Tel: (+55) 11 3732 3100  
Fax: (+55) 11 5522 9895

### Canada

Free Tel: 1800 565 1400  
Free Fax: 1800 265 3858  
Tel: (+1) 905 829 9500  
Fax: (+1) 905 829 9292

### Chile

Tel: (+56) 2 495 7395  
Fax: (+56) 2 495 7396

### People's Republic of China

Free Tel: 800 819 3336  
Tel: (+86) 21 6141 5566  
Fax: (+86) 21 6141 5567

### Czech Republic

Tel: (+420) 246 003 200  
Fax: (+420) 246 003 291

### Denmark

Tel: (+45) 43 56 59 00  
Fax: (+45) 43 56 59 05

### Finland

Tel: (+358) 9 350 9250  
Fax: (+358) 9 350 92555

### France

Free Tel: 0800 211 408  
Free Fax: 0800 031 052  
Tel: (+33) 474 82 28 88  
Fax: (+33) 474 95 68 08

### Germany

Free Tel: 0800 51 55 000  
Free Fax: 0800 64 90 000  
Tel: (+49) 89 6513 0  
Fax: (+49) 89 6513 1169

### Hungary

Tel: (+36) 1 235 9055  
Fax: (+36) 1 235 9068

### India

#### Telephone

Bangalore: (+91) 80 6621 9400  
New Delhi: (+91) 11 4358 8000  
Mumbai: (+91) 22 4087 2364  
Pune: (+91) 20 4146 4700  
Hyderabad: (+91) 40 3067 7450  
Kolkata: (+91) 33 4013 8000

#### Fax

Bangalore: (+91) 80 6621 9550  
New Delhi: (+91) 11 4358 8001  
Mumbai: (+91) 22 2579 7589  
Pune: (+91) 20 4146 4777  
Hyderabad: (+91) 40 3067 7451  
Kolkata: (+91) 33 4013 8016

### Ireland

Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: +353 (0) 402 20370  
Fax: +353 (0) 402 20375

### Israel

Free Tel: 1 800 70 2222  
Tel: (+972) 8 948 4222  
Fax: (+972) 8 948 4200

### Italy

Free Tel: 800 827 018  
Tel: (+39) 02 3341 7310  
Fax: (+39) 02 3801 0737

### Japan

Tel: (+81) 3 5796 7300  
Fax: (+81) 3 5796 7315

### Korea

Free Tel: (+82) 80 023 7111  
Free Fax: (+82) 80 023 8111  
Tel: (+82) 31 329 9000  
Fax: (+82) 31 329 9090

### Luxembourg

Tel: (+32) 3 899 1301  
Fax: (+32) 3 899 1311

### Malaysia

Tel: (+60) 3 5635 3321  
Fax: (+60) 3 5635 4116

### Mexico

Free Tel: 01 800 007 5300  
Free Fax: 01 800 712 9920  
Tel: (+52) 722 276 1600  
Fax: (+52) 722 276 1601

### The Netherlands

Tel: (+31) 78 620 5411  
Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Norway

Tel: (+47) 23 17 60 00  
Fax: (+47) 23 17 60 10

### Poland

Tel: (+48) 61 829 01 00  
Fax: (+48) 61 829 01 20

### Portugal

Free Tel: 800 202 180  
Free Fax: 800 202 178  
Tel: (+351) 21 924 2555  
Fax: (+351) 21 924 2610

### Russia

Tel: (+7) 495 621 5828  
Fax: (+7) 495 621 6037

### Singapore

Tel: (+65) 6779 1200  
Fax: (+65) 6779 1822

### Slovakia

Tel: (+421) 255 571 562  
Fax: (+421) 255 571 564

### South Africa

Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: (+27) 11 979 1188  
Fax: (+27) 11 979 1119

### Spain

Free Tel: 900 101 376  
Free Fax: 900 102 028  
Tel: (+34) 91 661 99 77  
Fax: (+34) 91 661 96 42

### Sweden

Tel: (+46) 8 742 4200  
Fax: (+46) 8 742 4243

### Switzerland

Free Tel: 0800 80 00 80  
Free Fax: 0800 80 00 81  
Tel: (+41) 81 755 2511  
Fax: (+41) 81 756 5449

### Thailand

Tel: (+66) 2 126 8141  
Fax: (+66) 2 126 8080

### United Kingdom

Free Tel: 0800 717 181  
Free Fax: 0800 378 785  
Tel: (+44) 01747 833 000  
Fax: (+44) 01747 833 574

### United States

Toll-Free: 800 325 3010  
Toll-Free Fax: 800 325 5052  
Tel: (+1) 314 771 5765  
Fax: (+1) 314 771 5757

### Vietnam

Tel: (+84) 8 3516 2810  
Fax: (+84) 8 6258 4238

### Internet

[sigma-aldrich.com](http://sigma-aldrich.com)

*Enabling Science to  
Improve the Quality of Life*

Order/Customer Service: [sigma-aldrich.com/order](http://sigma-aldrich.com/order)  
Technical Service: [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
Development/Custom Manufacturing Inquiries **S AFC**® [safcglobal@sial.com](mailto:safcglobal@sial.com)  
Safety-related Information: [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

World Headquarters  
3050 Spruce St.  
St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)