


# Microbiology Focus

Volume 6.2, 2014

**Fluka**<sup>®</sup>  
Analytical

## Microbiology Proficiency Testing



*How correct are your results? ISO 17025 accredited laboratories must regularly undergo proficiency testing (PT) to demonstrate their competency. Correct qualitative and quantitative results in a PT assure the labs that they are doing things the right way.*

Contents	
Microbiology	
Proficiency Testing	..... 2
New Photo Competition	.... 6
Confirmation of Bacteria by Enzymes	..... 6

# Microbiology proficiency testing

Jvo Siegrist, Product Manager Microbiology — [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)

We are living in a world with an increasing amount of regulation. However, in the field of microbiology, there is still some uncertainty with regard to quality control parameters.

Microbial safety is a major issue for water, plants, and the food and beverage industry, but there is also a great deal of pressure to reduce costs while increasing productivity. Currently, microbiology provides a broad range of recommendations to satisfy the regulations. Even though more methods are known and available on the market, the potential and demand for more rapid methods with greater accuracy is very high. On the one hand, the general trend is toward more standardization, but on the other hand, we have non-regulated new steps, new situations, new samples and new problems. Organizations such as ISO, AFNOR, UKAS, ASTM and ILAC propose guidelines and offer support regarding standardization in microbiology quality control.

The problems in a microbiology QC lab are broad, and compared to chemical analysis, there are more unanswered questions with regard to unknown and uncontrolled variables. It is sometimes difficult to decide if results obtained are correct or if they are the result of errors or a natural phenomenon.

## Here are some examples of problems:

- Large deviations in analysis
- Discrepancies between labs
- Difficult to compare results because of many parameters, different methods and different reporting of the results
- Confusion in determining the most suitable methods
- Validation takes a lot of time
- Human error, because of handling, calculation and reporting
- Failures tracing back to equipment, culture media, test, etc.

With the current trend toward standardization, microbiology is moving to the next level, especially due to the fact that knowledge about microbiology quality control has increased over the last 20 years. The methods employed are more accurate and labs do a better job of self-regulation.

## The following, however, needs to occur to ensure even more reliable results in microbiology quality control:

- Standardization of methods by following ISO, UKAS, EU Regulation, FDA or other National Accreditation bodies
- Good, defined and stable performance of tests in reference to qualitative and quantitative determinations
- Using standards like certified reference materials (reference strains)
- Participating in regular proficiency testing

## Did you know...

### More than 5% of false negatives pertain to the most important pathogens?

According to the American Proficiency Institute, a study with the four most common food pathogens - *E.coli* O157:H7, *Salmonella* spp., *Listeria monocytogenes*, and *Campylobacter* spp., showed that the average percentage of false negative results was consistently above 5.0% for all four pathogens throughout the study period (data collection over 11 years). Mainly, the issues are seen with *Campylobacter* spp.



Figure 1: Culture plates.

One basic but vital tool needed to improve the accuracy of testing is the use of **standards**. In the UKAS LAB 31 2nd edition, it is stated that control strains should be ISO 17025 certified if available. They should be comprised of material from a recognized national culture collection or a reference materials producer like Sigma-Aldrich, which produces under ISO guide 34 with strains from ATCC and NCTC.

In the case of a microbiology lab, standard means the need for both microorganism standards or a control strain, as well as standards for calibration of equipment, such as balances or incubators. Reference cultures are not only required for testing the performance of traditional culture methods and new techniques, but also for validation of new methods and to confirm the competency of the lab. It is also possible to use reference strains, which are derivatives of national or international reference cultures, as long as it can be proven that the relevant properties for the application still exist. According to ISO 11133-1, it is possible to culture reference strains for a single passage to produce reference stock cultures which are then controlled for purity and for biochemical tests. They should be stored in a freezer or in a freeze-dried form in small aliquots; however, defrosted cultures should not be refrozen. It is preferable for working strain cultures to be made out of stock cultures, and they should not be subcultured again.




Figure 2: PT samples.

### Proficiency Testing (PT) and ISO 17025

PT is an effective tool to help laboratories assure themselves and the accreditation bodies that the results that they report are correct, and to verify the effectiveness of the accreditation process. It is a comparison of results between several labs and demonstrates the technical competence of laboratories. ISO/IEC 17025 describes the general requirements for the competence of testing and calibration laboratories. It is specific as to the requirements for competence, and it applies directly to those organizations that produce testing and calibration results.

**Principle:** The labs getting samples of known, but undisclosed content, go through the routine procedures. As a result, the testing laboratory gets an assurance of their performance by an independent, external assessment.

**The philosophy of ISO/IEC 17025:** The same sample at different times, from different analysis, and from many different laboratories should reflect an agreeable result. Participating in a robust proficiency testing scheme not only gives laboratory managers confidence in their laboratory equipment, methodologies and laboratory staff, but also provides assurance that the laboratory is delivering the quality of results demanded by its customers. The accreditation according to ISO 17025 is a certification that guarantees calibrations and PT schemes are regularly performed – a true seal of approval!

### How does a PT process work?

In the flow chart of **Figure 4**, the process of a PT cycle from Sigma-Aldrich is shown. It starts as a program and the PT organization assumes the coordination for the participating laboratories. First, the registration form has to be filled in as shown in **Figure 3**. For registration to a PT program go to our web site: [sigma-aldrich.com/proficiencytesting](http://sigma-aldrich.com/proficiencytesting). Then the kits are sent out to the labs. They contain different microorganisms in different concentrations (or blanks) with a sample matrix. The labs complete their testing and submit the results to the PT organization. The PT organization collects all results and performs a statistical analysis. A report is then generated and sent back to the participating laboratories.

The basis of these tests builds the Vitroids™ technology, which allow us to make a PT material with stable and exact colony forming units (CFU) and a narrow standard deviation (down to ± 4% on the level of 100 CFU). These microorganism standards are produced and sold as certified reference materials and are



Figure 3: Registration form.




Figure 4: PT Flow Chart.



Figure 5: Vitroids™ Discs.

used for PT schemes. They are made out of reliable reference strains from ATCC and NCTC, are produced under ISO guide 34, and the CFU value is certified under ISO 17025. The organisms are immobilized in a disc in a possible range of 30 to 10<sup>9</sup> CFU per disc.

The discs are easy to use since they can be placed directly in water, diluent, broth, or even on agar plates. The Vitroids™ contain highly viable bacteria and when placed in contact with media, they dissolve rapidly and start to grow without a lag-phase. The viability of the CFU in a disc is stable for at least one

year (for most organisms, more than two years) when kept under refrigeration (-20 °C). It is also acceptable for the product to be briefly transported at ambient temperature. Each disc is packed individually with some desiccant and then sealed in mylar foil. Each package comes with a comprehensive certificate of analysis reporting the CFU and standard deviation.

Vitroids help microbiologist have reliable results, save a lot of time (laboratory, documentation), and lower their costs. For additional info visit: [sigma-aldrich.com/vitroids](http://sigma-aldrich.com/vitroids)

PT Cat. No.	QC Cat. No.	Description
MIC020		<i>Clostridium perfringens</i> in Water
MIC014	QCMIC014	<i>E. coli</i> - Sludge
MIC107	QCMIC107	<i>E. coli</i> in Drinking & Surface Water - Quantitative - WS 3 Levels
MIC207	QCMIC207	<i>E. coli</i> in Drinking & Surface Water - Quantitative- WS 4 Levels
MIC007	QCMIC007	<i>E. coli</i> in Drinking and Surface Water - Quantitative - WS
MIC022		<i>E. coli</i> in Seawater - Quantitative WP
MIC122		<i>E. coli</i> in Seawater - Quantitative WP 3 Levels
MIC003	QCMIC003	<i>E. coli</i> in Water - Quantitative WP
MIC103	QCMIC103	<i>E. coli</i> in Water - Quantitative WP 3 Levels
MIC009	QCMIC009	<i>E. coli</i> Quantitative - Soil
MIC109	QCMIC109	<i>E. coli</i> Quantitative - Soil 3 Levels
MIC209	QCMIC209	<i>E. coli</i> Quantitative - Soil 4 Levels
MIC015	QCMIC015	Fungi and Yeast - WS
MIC004	QCMIC004	<i>Legionella</i> in Water - WP
MIC019	QCMIC019	<i>Listeria</i> - WP
MIC008	QCMIC008	<i>Pseudomonas aeruginosa</i> - WS
MIC021		<i>Salmonella</i> - P/A 5 samples
MIC013	QCMIC013	<i>Salmonella</i> - Sludge
MIC113	QCMIC113	<i>Salmonella</i> - Sludge 3 Levels
MIC213	QCMIC213	<i>Salmonella</i> - Sludge 4 Levels
MIC106	QCMIC106	<i>Salmonella</i> - WS 3 Levels
MIC206	QCMIC206	<i>Salmonella</i> - WS 4 Levels
MIC006	QCMIC006	<i>Salmonella</i> for Drinking/Surface Water - WS
MIC012	QCMIC012	Standard Plate Count - WP
MIC112	QCMIC112	Standard Plate Count - WP 3 Levels
MIC212	QCMIC212	Standard Plate Count - WP 4 Levels
MIC002	QCMIC002	Standard Plate Count - WS
MIC102	QCMIC102	Standard Plate Count - WS 3 Levels
MIC202	QCMIC202	Standard Plate Count - WS 4 Levels
MIC011	QCMIC011	<i>Streptococcus/Enterococcus</i> - Drinking & Surface Water
MIC111	QCMIC111	<i>Streptococcus/Enterococcus</i> - Drinking & Surface Water 3 Levels
MIC211	QCMIC211	<i>Streptococcus/Enterococcus</i> - Drinking & Surface Water 4 Levels
MIC123		<i>Streptococcus/Enterococcus</i> - Seawater WP 3 Levels
MIC105	QCMIC105	<i>Streptococcus/Enterococcus</i> - WP 3 Levels
MIC023		<i>Streptococcus/Enterococcus</i> in Seawater - WP
MIC205	QCMIC205	Total & Fecal <i>Streptococcus/Enterococcus</i> - WP 4 Levels
MIC005	QCMIC005	Total and Fecal <i>Streptococcus/Enterococcus</i> - WP
MIC016	QCMIC016	WS-Enterococci -Sample (1-10)
MIC001	QCMIC001	WS-Microbiological - Sample (1-10)

**Table 1:** Proficiency testing schemes (PT) and Quality Check Sets (QC) for internal validation. The PTs and QCs are used to test the performance of microbial quality control in water supply (WS), water pollution (WP) and other matrices. The levels show how many concentrations (always in duplicate) are present in the set. More detailed information on each PT or QC set can be found on the web ([sigma-aldrich.com](http://sigma-aldrich.com)). Be aware that some products need a buffer for hydrolyzing; we recommend using Z699489 Phosphate Buffer, Magnesium Chloride, volume 72 x 99 mL.



Vitroids™	Origin	Strain No.	CFU	Cat. No.
<i>Aspergillus brasiliensis</i>	ATCC	16404™	80	RQC15003
<i>Bacillus subtilis</i>	ATCC	6633™	10'000	RQC02258
<i>Bacillus subtilis</i>	ATCC	6633™	80	RQC16003
<i>Candida albicans</i>	ATCC	10231™	80	RQC14003
<i>Candida albicans</i>	ATCC	10231™	1'000	RQC14007
<i>Candida albicans</i>	ATCC	10231™	10'000	RQC14008
<i>Citrobacter freundii</i>	ATCC	8090	200	RQC02105
<i>Clostridium perfringens</i>	NCTC	10240	30	RQC02351
<i>Clostridium perfringens</i>	NCTC	10240	200	RQC02355
<i>Clostridium perfringens</i>	NCTC	10240	500	RQC20106
<i>Clostridium sporogenes</i>	ATCC	19404™	80	RQC19003
<i>Enterobacter aerogenes</i>	ATCC	13048™	10'000	RQC01658
<i>Enterobacter aerogenes</i>	ATCC	13048™	100	RQC01654
<i>Enterobacter aerogenes</i>	ATCC	13048™	50	RQC01652
<i>Enterobacter aerogenes</i>	ATCC	13048™	200	RQC01655
<i>Enterobacter aerogenes</i>	ATCC	13048™	1'000	RQC01657
<i>Enterococcus cloacae</i>	ATCC	35030™	50	RQC21102
<i>Enterococcus faecalis</i>	ATCC	19433™	50	RQC01772
<i>Enterococcus faecalis</i>	ATCC	19433™	100	RQC01774
<i>Enterococcus faecalis</i>	ATCC	19433™	200	RQC01775
<i>Enterococcus faecalis</i>	ATCC	19433™	1'000	RQC01777
<i>Escherichia coli</i>	ATCC	11775™	1'000	RQC01707
<i>Escherichia coli</i>	ATCC	11775™	200	RQC01705
<i>Escherichia coli</i>	ATCC	11775™	50	RQC01702
<i>Escherichia coli</i>	ATCC	8739™	80	RQC11003
<i>Escherichia coli</i>	ATCC	11775™	10'000	RQC01708
Heterotrophic Organisms	—	—	100	RQC02504
<i>Fluoribacter bozemanae</i>	NCTC	11368	50'000	RQC02908
<i>Legionella pneumophila</i> (serogroup 1)	NCTC	12821	100'000	RQC02008
<i>Listeria monocytogenes</i>	ATCC	19115™	30	RQC01901
<i>Pseudomonas aeruginosa</i>	ATCC	9027™	30	RQC02202
<i>Pseudomonas aeruginosa</i>	ATCC	9027™	100	RQC02204
<i>Pseudomonas aeruginosa</i>	ATCC	9027™	50	RQC12002
<i>Pseudomonas aeruginosa</i>	ATCC	9027™	200	RQC12005
<i>Pseudomonas aeruginosa</i>	ATCC	9027™	1'000	RQC12007
<i>Salmonella enterica</i> subsp. Enterica serovar Typhimurium	ATCC	14028	50	RQC17002
<i>Salmonella enterica</i> subsp. Enterica serovar Abony	NCTC	6017	80	RQC18003
<i>Salmonella goldcoast</i>	NCTC	13175	30	RQC02301
<i>Staphylococcus aureus</i> susp. Aureus	ATCC	6538	1'000	RQC13007
<i>Staphylococcus aureus</i> susp. Aureus	ATCC	6538	200	RQC13005
<i>Staphylococcus aureus</i> susp. Aureus	ATCC	6538	50	RQC13002
Vitroids™ Blank	—	—	0	RQC0001

Table 2: Range of Vitroids

# Confirmation of bacteria by specific enzymes

Jvo Siegrist, Product Manager Microbiology — [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)

**A quick and easy confirmation of bacteria is always recommended at the end of a bacteria detection process.**

Confirmation steps are really important in microbiology, and there is a wide range of methods that can be employed. Often another differential medium or method is used, such as a biochemical or immunological test. In most cases, it is important that the test be performed within a short time in order to release the results as soon as possible, and it should not be too costly. Therefore, Sigma-Aldrich focused their effort on bacteria specific enzymes and a simple detection system with specific chromogenic and fluorogenic substrates. Test cards were developed, where the reaction is visible to the naked eye, and the result is obtained within 10 minutes. A nice product range has been provided to confirm *E. coli*, *Salmonella*, *Neisseria gonorrhoea*, *Staphylococcus aureus*, *Enterococci*, Group A *Streptococci*, *Total coliform*, *Fecal coliform* and *Gram negatives*.

These tests utilize a specific substrate which, when hydrolyzed by a specific enzyme of the target organism (during peptide hydrolysis), produces a blue/white fluorescence, or purple/blue color upon the addition of Reagent B (color developer) when applicable.

## Content of Kit

The test set contains 12 test cards with 4 test spots on each card (sufficient for 48 tests) and a reagent A (buffer), and in some tests, a reagent B with color developer is provided. Each kit contains a product insert. All materials can be stored at room temperature (away from direct sunlight). Refrigeration storage does not harm the test. The shelf life of the test is more than 2 years. Cards which are no longer white should not be used anymore and could make the results invalid.

## Validation

Over 500 isolates were tested with these test cards, with a correlation of greater than 99% when compared to traditional biochemical testing.

## Quality Control/Incubation

When using the test, it should be checked each time with positive and negative controls, using known stock strains of the target organism and another strain of bacteria that is not the target organism. Ideally, cultured bacteria should be no more than 18 hours old and grown on an agar plate specified in **Table 1**, otherwise the target enzyme may no longer be present (it is only produced during the exponential growth phase). Using other media

## 5th Microbiology photo competition

**Win an android™ tablet in the fluka microbiology competition.**

This photography competition is sponsored by Sigma-Aldrich with the aim of encouraging microbiologists to promote something about their work and their science. The best photographic entries will win prizes such as an Android™ tablet PC, a Giant Microbe, a Swiss army knife and a USB stick. The winning images will be published in Microbiology Focus, and the best one will be featured on the cover.

### Rules of the Competition and Conditions of Entry

1. The competition is open to all residents worldwide.
2. Entries should illustrate any microorganisms (living or dead) or a microbiologist in action at work
3. Picture size should be at least 400 dpi and 90 x 120 mm (max 5 MB). The file format must be in jpg, tiff or pdf!
4. The entries will be judged on:
  - Clarity of presentation
  - Composition
  - Illumination and contrast
  - Congruency of subject matter and title of photograph
  - Scientific interest and relevance
  - Originality



5. Winning entries will be retained by Sigma-Aldrich, who will have sole rights of publication, reproduction and display.
6. Closing date will be 31st August 2014
7. Entries after the closing date will not be considered. Entries that are incomplete, illegible, mutilated, altered or not complying exactly with the instructions and theme may be disqualified.
8. Decisions of the judges in all matters affecting the competition will be final and legally binding.

### The competition will be judged by:

**Dr. Lars Fieseler**, Zurich University of Applied Sciences - ZHAW, Supervisor, Department Microbiology

**Prof. Mohammad Manafi**, Medical University of Vienna, Head of the Department for Food Hygiene

**Jvo Siegrist**, Sigma-Aldrich, Product Manager, Microbiology

### Method of Entry

There is no entry fee, but for each entry, an entry form must be completed (three entries at the most).


Entry form available from

**[sigma-aldrich/mibi-competition](http://sigma-aldrich/mibi-competition)**


could also influence the metabolism and cause false results, so it is recommended that only the specified media in **Table 1** be used.

### Procedure

Use a culture from an agar plate of the types specified in Table 1. Other types of growth media may not allow the target organism to produce the target enzyme to be detected by the test card. Ideally, the growth should be between 14 to 18 hours old before inoculating the test card. After 24 hours of growth, most bacteria go into a stationary phase where reduced levels of enzymes are present, which may affect the results. If a test of a negative result is done after 24 hours of growth and a concern remains, repeat the culture and test within 14 to 18 hours. (Note: The incubation temperature should be 44 °C for fecal coliform, and 38.5 °C for all other types of bacteria.)



**Figure 1:** Blue/white fluorescence on an Enzyme Confirmation Card.



**Figure 2:** Color reaction on an Enzyme Confirmation Card.

1. Select one test spot as your Control Spot, and other test spots as your Sample Spot or positive control. Add 1 drop of Reagent A to the filter membrane area of the Control Spot and to the filter membrane area of each Sample Spot.
2. Select colonies that morphologically resemble the target organism from the first growth plate and touch the tops of 1-2 colonies with a loop, inoculating needle, swab, or wooden applicator. Smear the colonies onto the filter membrane area of the Sample Spot ONLY.

Set aside the inoculated card at room temperature for 5-10 minutes.

Part Number	Target Organism	Inspection			Suitable Growth Agar Media										
		Reagent B	Method	Shown Color if Positive	B	M	E	Sal	S.Sal	GN	EC	110	LT	CH	
*75444	<i>E. coli</i>	No	UV Light	Blue/White Fluorescence	X	X					X	X			
*77643	Total Coliform	No	UV Light	Blue/White Fluorescence	X	X					X			X	
*40926	Fecal Coliform	No	UV Light	Blue/White Fluorescence	X	X					X			X	
*55283	<i>Salmonella</i>	No	UV Light	Yellow Fluorescence	X			X	X						
*56305	<i>Enterococcus</i>	Yes	Visual	Purple	X		X								
*74203	Gram+/Gram-Differentiation	Yes	Visual	Purple (indicate Gram-)	X						X				
*92598	<i>Neisseria gonorrhoea</i>	Yes	Visual	Blue											X
*77701	Group A <i>Streptococcus</i>	Yes	Visual	Purple	X										
*80031	<i>Staphylococcus aureus</i>	Yes	Visual	Purple	X								X		

\*Currently not available in the USA.

**Table 1:** Range of Enzyme Confirmation Test Cards/Soon also available in USA

#### Index List of Media:


B = Blood Agar (e.g., Fluka 70133)  
M = MacConkey Agar (e.g., Fluka 70143)  
E = *Enterococcus* Selective Agar (e.g., Fluka 45183)  
Sal = *Salmonella* Agar (W/brilliant green) (e.g., Brilliant Green Agar, modified Fluka 70134?)  
S.Sal = Selective *Salmonella* Agar (e.g. Fluka 05538 or 738419?)

## Did you know...

### Bacteria do not waste time and energy?

Generally, bacteria do not synthesize enzymes unless the substrates are present and they are needed. The regulation of the enzyme production is an interplay of substrates, inducers and inhibitors.

**Figure 3:** ATP synthase is an enzyme that provides energy for cells by synthesizing adenosine triphosphate.



**Figure 3**

A) UV detection method (see **Table 1**)

3. Look at the test card under a long wave (~360nm) UV light. A positive test will be indicated by a fluorescence light (see **Table 1**) around the smeared colonies in the Sample Spot over a blue background. The test is negative if no fluorescence colonies appear in the Sample Spot (which will look like the Control Spot). A positive test is indicative of the target organism. The fluorescence of a positive result should remain on the Sample Spot for more than 24 hours.

B) Color detection method (see **Table 1**)

4. Add 1 drop of Reagent B to the filter membrane area of the Sample Spot and to the filter membrane area of each Sample Spot.
5. A purple or blue (as stated in **Table 1**) color will immediately form on and around the deposited colonies on the Sample Spot in the presence of the hydrolyzed substrate, indicating a Positive test of the target organism. The test is negative if no purple or blue color appears in the Sample Spot. The test result must be read within 1 minute after the addition of Reagent B, otherwise the test result may be invalid. If purple or blue color shows on the Control Spot, then re-test.

### Limitations

Occasionally, some species of other microorganisms may produce small amounts of enzyme, which can produce a positive test from culture. This is extremely rare. The test detects bacterial enzymes, which may dissolve in liquids. Serial dilutions of bacteria may not give a positive reaction.

## Sigma-Aldrich® Worldwide Offices

### Argentina

Free Tel: 0810 888 7446  
Tel: (+54) 11 4556 1472  
Fax: (+54) 11 4552 1698

### Australia

Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Austria

Tel: (+43) 1 605 81 10  
Fax: (+43) 1 605 81 20

### Belgium

Tel: (+32) 3 899 13 01  
Fax: (+32) 3 899 13 11

### Brazil

Free Tel: 0800 701 7425  
Tel: (+55) 11 3732 3100  
Fax: (+55) 11 5522 9895

### Canada

Free Tel: 1800 565 1400  
Free Fax: 1800 265 3858  
Tel: (+1) 905 829 9500  
Fax: (+1) 905 829 9292

### Chile

Tel: (+56) 2 495 7395  
Fax: (+56) 2 495 7396

### People's Republic of China

Free Tel: 800 819 3336  
Tel: (+86) 21 6141 5566  
Fax: (+86) 21 6141 5567

### Czech Republic

Tel: (+420) 246 003 200  
Fax: (+420) 246 003 291

### Denmark

Tel: (+45) 43 56 59 00  
Fax: (+45) 43 56 59 05

### Finland

Tel: (+358) 9 350 9250  
Fax: (+358) 9 350 92555

### France

Free Tel: 0800 211 408  
Free Fax: 0800 031 052  
Tel: (+33) 474 82 28 88  
Fax: (+33) 474 95 68 08

### Germany

Free Tel: 0800 51 55 000  
Free Fax: 0800 64 90 000  
Tel: (+49) 89 6513 0  
Fax: (+49) 89 6513 1169

### Hungary

Tel: (+36) 1 235 9055  
Fax: (+36) 1 235 9068

### India

#### Telephone

Bangalore: (+91) 80 6621 9400  
New Delhi: (+91) 11 4358 8000  
Mumbai: (+91) 22 4087 2364  
Pune: (+91) 20 4146 4700  
Hyderabad: (+91) 40 3067 7450  
Kolkata: (+91) 33 4013 8000

#### Fax

Bangalore: (+91) 80 6621 9550  
New Delhi: (+91) 11 4358 8001  
Mumbai: (+91) 22 2579 7589  
Pune: (+91) 20 4146 4777  
Hyderabad: (+91) 40 3067 7451  
Kolkata: (+91) 33 4013 8016

### Ireland

Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: +353 (0) 402 20370  
Fax: +353 (0) 402 20375

### Israel

Free Tel: 1 800 70 2222  
Tel: (+972) 8 948 4222  
Fax: (+972) 8 948 4200

### Italy

Free Tel: 800 827 018  
Tel: (+39) 02 3341 7310  
Fax: (+39) 02 3801 0737

### Japan

Tel: (+81) 3 5796 7300  
Fax: (+81) 3 5796 7315

### Korea

Free Tel: (+82) 80 023 7111  
Free Fax: (+82) 80 023 8111  
Tel: (+82) 31 329 9000  
Fax: (+82) 31 329 9090

### Luxembourg

Tel: (+32) 3 899 1301  
Fax: (+32) 3 899 1311

### Malaysia

Tel: (+60) 3 5635 3321  
Fax: (+60) 3 5635 4116

### Mexico

Free Tel: 01 800 007 5300  
Free Fax: 01 800 712 9920  
Tel: (+52) 722 276 1600  
Fax: (+52) 722 276 1601

### The Netherlands

Tel: (+31) 78 620 5411  
Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Norway

Tel: (+47) 23 17 60 00  
Fax: (+47) 23 17 60 10

### Poland

Tel: (+48) 61 829 01 00  
Fax: (+48) 61 829 01 20

### Portugal

Free Tel: 800 202 180  
Free Fax: 800 202 178  
Tel: (+351) 21 924 2555  
Fax: (+351) 21 924 2610

### Russia

Tel: (+7) 495 621 5828  
Fax: (+7) 495 621 6037

### Singapore

Tel: (+65) 6779 1200  
Fax: (+65) 6779 1822

### Slovakia

Tel: (+421) 255 571 562  
Fax: (+421) 255 571 564

### South Africa

Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: (+27) 11 979 1188  
Fax: (+27) 11 979 1119

### Spain

Free Tel: 900 101 376  
Free Fax: 900 102 028  
Tel: (+34) 91 661 99 77  
Fax: (+34) 91 661 96 42

### Sweden

Tel: (+46) 8 742 4200  
Fax: (+46) 8 742 4243

### Switzerland

Free Tel: 0800 80 00 80  
Free Fax: 0800 80 00 81  
Tel: (+41) 81 755 2511  
Fax: (+41) 81 756 5449

### Thailand

Tel: (+66) 2 126 8141  
Fax: (+66) 2 126 8080

### United Kingdom

Free Tel: 0800 717 181  
Free Fax: 0800 378 785  
Tel: (+44) 01747 833 000  
Fax: (+44) 01747 833 574

### United States

Toll-Free: 800 325 3010  
Toll-Free Fax: 800 325 5052  
Tel: (+1) 314 771 5765  
Fax: (+1) 314 771 5757

### Vietnam

Tel: (+84) 8 3516 2810  
Fax: (+84) 8 6258 4238

### Internet

[sigma-aldrich.com](http://sigma-aldrich.com)

*Enabling Science to  
Improve the Quality of Life*

Order/Customer Service: [sigma-aldrich.com/order](http://sigma-aldrich.com/order)  
Technical Service: [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
Development/Custom Manufacturing Inquiries: [SAFC® safcglobal@sial.com](mailto:safcglobal@sial.com)  
Safety-related Information: [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

World Headquarters  
3050 Spruce St.  
St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)